

# EQUASHIELD® Pro Overview

Equashield® is a developer of Closed System Transfer Devices (CSTDs) for the safe handling of hazardous drugs to protect healthcare workers. Equashield's mechanically closed CSTD design relies on its unique pressure equalization behind the syringe and with its use of tight seal double membranes to keep connectors residual free despite multiple entries. Equashield® Pro is built around the Equashield CSTD and relies on its innovative design to redefine the automation workflow and optimize throughput.

Unlike traditional robot's reliance on robotic arms that mimic human motion with use of needle and syringe leading to potential for cytotoxic drips, spills, aerosols and vapor release, Equashield® Pro uses multi-station compounding approach to minimize both the motion and time it takes to prepare a dose much like a factory manufacturing line. Equashield® Pro is also capable of processing a broad variety of finished product dosage forms and access all size drug vials to be a truly versatile machine. Equashield® Pro is designed to make adoption of robotics in the IV room simple and affordable regardless of hospital size. It is easy to fit into hospital pharmacies and comparable in size to a standard Biological Safety Cabinet.



**Hazardous Drug Protection**  
with use of Equashield's leading Closed System Transfer Device tested for Leak Free connection and Air Tight closeness



**High Throughput** to meet your compounding needs during peak time, enabling you to optimize patient care without sacrificing scarce resources



**Medication Error Protection**  
with use of Equashield's proprietary dose verification system to prevent wrong drug, wrong dose and labeling errors during compounding



**Small Footprint** to fit into existing pharmacy space. Equashield Pro is comparable to the size of a standard Biological Safety Cabinet



**Highly Versatile** to produce a wide range of patient specific dose with syringes, IV bags and other finished dose forms



**Cost Effective Solution** to fit into any pharmacy budget. Additional cost avoidance from expensive renovation costs and expensive maintenance

## EQUASHIELD Pro Key Metrics

67

Large Vial  
Storage Capacity

8

Simultaneous  
Work Stations

65

Syringe Unit  
Storage Capacity

8

Shakers with 2  
Reconstitution  
Stations

60

Doses Compounded  
Per Hour

# EQUASHIELD®

Safe. Simple. Closed.



/equashield

# AUTOMATED COMPOUNDING of HAZARDOUS DRUGS

## EQUASHIELD® Pro



EQP-1601.1 ©2017 Equashield, LLC

EQUASHIELD, LLC | 99 Seaview Blvd. Port Washington, NY 11050  
855-EQUASHIELD | info@equashield.com | www.equashield.com


# A Progressive Choice in Automated Hazardous Drug Compounding

- ✓ Protection from Hazardous Drugs via the use of Clinically Proven CSTD System
- ✓ Highest Throughput in the Industry when preparing Patient Specific Doses by leveraging Equashield CSTD System
- ✓ Protection from Medication Errors during Compounding via the use of dose verification system



## Robotic Hazardous Drug Compounding With Built-In Safety of Industry Leading Closed System Transfer Device (CSTD)


# A Look Inside EQUASHIELD® Pro




# Small Footprint – Large Capacity

PRODUCT SPECIFICATION	PRODUCT REQUIREMENT
Size of Pro (WxDxH).....	2.223 x 1.353 x 2.261 (m), 87 x 54 x 89 (in)
Minimum Clearance (in).....	36 (front), 12 (Right), 12 (Left), 12 (Back)
Minimum Distance from Ceiling.....	0.15 (m), 6 (in)
Inner workspace of Equashield Pro (w x d).....	72 x 30 (in)
Weight (Empty).....	1,100 (Kg), 2425 (lb)
Floor Rating.....	>113 lb/ft2 (550 Kg/m2)
Doorway Clearance Necessary for Installation.....	0.978 (m), 38.5 (in)
Corridor Width Requirement (to mitigate 90 degree turns).....	1.485 (m), 58.5 (in)
Air Outflow Max.....	3100 M3/H
Power Supply.....	190V, 20A, 3P
Networking Interface.....	2x RJ45 (LAN) with HL7 Interface

Front View



Side View



## PRODUCT CAPACITY

**VIALS**  
Storage Capacity: Over 67 large vials (100ml),  
Sizes: 1ml to 100ml,  
Manufacturer: Most commercially available vials of 13, 20 and 28mm neck size

**SYRINGES**  
Conveyor Capacity: 65, Sizes: 1ml to 60ml,  
Manufacturer: Equashield Syringe Units

**IV BAGS**  
Simultaneous Workstations: 8, Sizes: 50ml to 1000 ml,  
Manufacturer: Any commercially available IV Bags with ISO compatible IV Port

Top View

