

AORFIX™

Endovascular Stent Graft

**THE ONLY ONE APPROVED BY
FDA IN NECK ANGLES UP TO 90°**

AORFIX™ BODY


PROXIMAL DIAMETER (mm)

D1: 24 / 27 / 31

D1

GRAFT BODY
LENGTH (mm)
L1: 81 / 96 / 111 / 126

L1



CONTRALATERAL LEG

LENGTH (mm)

L3: L1 > 81


56 / 64 / 73 / 81 / 90 / 98 / 106

12mm

12mm


OVERLAP
40mm

L3



PROXIMAL EXTENDER

D5



L5

DIAMETER (mm)

D5: 24 / 27 / 31

LENGTH (mm)

L5: 38

DISTAL EXTENDER

D6



L6

DIAMETER (mm)


D6: 10 / 12 / 14 / 16 / 18 / 20

LENGTH (mm)

L6: 51 / 82

CONVERTER

D4



L4

DIAMETER (mm)

D4: 25 / 27 / 29 / 31

LENGTH (mm)

L4: 72